

## Ostomy Accessories

### Hollister

#### Medical Adhesive

- Medical-grade silicone adhesive is unaffected by moisture
- Can be used on sensitive skin and allows skin to breathe
- Remains flexible and comfortable and provides a secure and reliable seal
- Easy-to-use spray allows application directly to appliance prior to application to skin



M

CIN	Mfr. No.	Description	Size	Qty.
3278199	7730	Adhesive	3.2 oz.	1/ea

HCPC Code: A4364

### ConvaTec®

#### Stomahesive® Paste

- Hydrocolloid-based protective skin barrier
- Use as a filler in uneven skin surfaces to help increase ostomy system wear-time and protect skin



M

CIN	Mfr. No.	Description	Size	Qty.
1019124	183910	Paste	2 oz.	1/ea

HCPC Code: A4406

### ConvaTec®

#### Stomahesive® Protective Powder

- Helps form a protective barrier for excoriated or weeping skin



M

CIN	Mfr. No.	Description	Size	Qty.
1022003	25510	Powder	1 oz.	1/ea

HCPC Code: A4371

### Hollister

#### Adapt® Paste

- Used to caulk and fill skin crevices to improve the fit of a skin barrier
- Low-alcohol paste helps to minimize stinging



M

CIN	Mfr. No.	Description	Size	Qty.
3406469	79300	Paste	2 oz.	1/ea

HCPC Code: A4406

### Hollister

#### Adapt® Stoma Powder

- Protects raw and weeping skin by creating a sticky gel
- Use sparingly to absorb moisture or exudate from skin prior to placing a skin barrier on peristomal skin
- Convenient to use; squeeze bottle allows aeration of powder onto skin; bottle has window for visualizing contents



M

CIN	Mfr. No.	Description	Size	Qty.
2982098	7906	Powder	1 oz.	1/ea

HCPC Code: A4371

### Hollister

#### Adapt® Convex Barrier Ring

- Can be cut, bent and stacked together to improve the fit of the skin barrier
- Used to adjust skin barrier thickness for deeper convexity
- Used to create a custom fit, including oval-shaped convexity



M

CIN	Mfr. No.	Description	Size	Qty.
3406832	79530	Barrier ring	1 3/8 in. inner diameter	10/ct

HCPC Code: A4385

### Hollister

#### Adapt® Flat Barrier Ring

- Sting-free alternative to paste
- Can be stretched and molded to create custom shapes
- Can be cut, bent and stacked together to improve the fit of the skin barrier
- For individuals with sensitive skin or limited dexterity
- Prolongs skin barrier wear time when used under a pouch or two-piece skin barrier



M

CIN	Mfr. No.	Description	Size (outer)	Qty.
3088036	7805	Barrier ring	2 in.	10/ct
3085206	7806	Barrier ring	4 in.	10/ct

HCPC Code: A4385

M May be reimbursable by Medicare Part B

Available through Cardinal Health Drop Ship program

Note: Products may not be available in all Cardinal Health Distribution Centers

## Ostomy Accessories

ConvaTec®

### Eakin Cohesive® Seals

- Can be molded to a variety of shapes to help prevent leakage around the stoma by acting as filler for uneven skin surfaces



CIN	Mfr. No.	Description	Size	Qty.
2572642	839002	Small	2 in.	20/ct
2572634	839001	Large	4 in.	10/ct

HCPC Code: A4385

ConvaTec®

### Stomahesive® Skin Barrier

- A multi-purpose skin barrier designed to help protect peristomal skin
- Gentle on skin yet holds firmly to skin surface



CIN	Mfr. No.	Description	Size	Qty.
1191105	21712	Nonsterile, no starter hole	4 x 4 in.	5/ct

HCPC Code: A4362

ConvaTec®

### AllKare® Adhesive Remover Wipes

- Gently removes all adhesives including skin barriers, tapes and hydrocolloid dressings



CIN	Mfr. No.	Description	Qty.
2827244	37436	Wipes	50/ct
1678804	37443	Wipes	100/ct

HCPC Code: A4365

ConvaTec®

### AllKare® Protective Barrier Wipes

- Provides a barrier film layer on skin used under tapes, skin barriers, adhesives and hydrocolloid dressings
- Helps protect against irritation, excoriation and adhesive build-up



CIN	Mfr. No.	Description	Qty.
2827228	37439	Wipes	50/ct
1678788	37444	Wipes	100/ct

HCPC Code: A5120

Hollister

### Universal Remover Wipes for Adhesives and Barriers

- Gentle solvent removes adhesives and cements
- Effective for most adhesives and safe to use on adult skin
- Easy-to-use in a convenient wipe package



CIN	Mfr. No.	Description	Qty.
2881183	7760	Remover wipes	50/ct

HCPC Code: A4456

Hollister

### Adapt® Lubricating Deodorant

- Enables a new level of confidence and security with a proven odor-neutralizer and pouch lubricant that helps keep the pouch interior clean
- Lubrication eases emptying and helps prevent pouch static and sticking
- Deodorant works by neutralizing the odor; not a perfume or masking agent



CIN	Mfr. No.	Description	Size	Qty.
3593571	78500	Deodorant	8 oz.	1/ea

HCPC Code: A4394

Hollister

### m9™ Odor Eliminator Drops

- Highly effective in eliminating odors; safe to use as demonstrated over many years
- Convenient to use; place drops of m9™ into pouch before applying it
- Economical; a few drops per wearing is adequate to eliminate odor



CIN	Mfr. No.	Description	Size	Qty.
3157831	7715	Bottle	1 oz.	1/ea
3278181	7717	Bottle	8 oz.	1/ea

HCPC Code: A4394

ConvaTec®

### Ostomy Appliance Belt

- Provides an added sense of security
- Adjustable



CIN	Mfr. No.	Description	Size	Qty.
1019157	175507	Appliance belt	Up to 42 in.	1/ea

ConvaTec®

**ActiveLife® One-Piece Drainable Pouch Cut-to-Fit Stomahesive® with Tape Collar**

- A lightweight, flexible, one-piece system
- Combines the skin barrier and the pouch into one simple, easy-to-manage system
- Pouch allows for draining
- One-sided comfort panel
- Can accommodate ¾–2½ in.



M

CIN	Mfr. No.	Color	Barrier Opening	Size	Qty.
1039015	22771	Transparent	Cut-to-fit	12 in.	10/ct

HCPC Code: A5061

ConvaTec®

**ActiveLife® One-Piece Drainable Pouch Pre-Cut Stomahesive® with Tape Collar**

- A lightweight, flexible, one-piece system
- Combines the skin barrier and the pouch into one simple, easy-to-manage system
- Pouch allows for draining



M

CIN	Mfr. No.	Barrier Opening	Size	Qty.
<b>1-sided comfort panel, opaque</b>				
1038637	22757	¾ in.	12 in.	10/ct
1038652	22758	1 in.	12 in.	10/ct
1038686	22759	1¼ in.	12 in.	10/ct
1038777	22760	1½ in.	12 in.	10/ct
1038793	22761	1¾ in.	12 in.	10/ct
<b>1-sided comfort panel, transparent</b>				
1038884	22765	1 in.	12 in.	10/ct
1038900	22766	1¼ in.	12 in.	10/ct
1038934	22767	1½ in.	12 in.	10/ct
1038967	22768	1¾ in.	12 in.	10/ct
1038975	22769	2 in.	12 in.	10/ct
<b>1-sided comfort panel, opaque</b>				
1046630	22750	¾ in.	10 in.	10/ct
1046648	22751	1 in.	10 in.	10/ct
1038561	22752	1¼ in.	10 in.	10/ct
1038579	22753	1½ in.	10 in.	10/ct

HCPC Code: A5061

ConvaTec®

**ActiveLife® One-Piece Convex Drainable Pouch Pre-Cut Durahesive® with Tape Collar**

- A lightweight, flexible, one-piece system
- Combines the skin barrier and the pouch into one simple, easy-to-manage system
- Convex barrier is good for people with a recessed or flat stoma
- Pouch allows for draining



M

CIN	Mfr. No.	Barrier Opening	Size	Qty.
<b>1-sided comfort panel, transparent</b>				
1129873	175777	¾ in.	12 in.	5/ct
1129964	175778	⅞ in.	12 in.	5/ct
1130129	175779	1 in.	12 in.	5/ct
1130152	175780	1⅛ in.	12 in.	5/ct
1130178	175781	1¼ in.	12 in.	5/ct
1130194	175782	1⅜ in.	12 in.	5/ct
1130244	175783	1½ in.	12 in.	5/ct
3581410	125352	¾ in.	12 in.	10/ct
3581428	125353	⅞ in.	12 in.	10/ct
3581444	125354	1 in.	12 in.	10/ct
3581451	125355	1⅛ in.	12 in.	10/ct
3581469	125356	1¼ in.	12 in.	10/ct

HCPC Code: A4390

ConvaTec®

**ActiveLife® Convex One-Piece Urostomy Pouch Pre-Cut Durahesive®**

- A lightweight, flexible, one-piece system
- Combines the skin barrier and the pouch into one simple, easy-to-manage system
- Convex barrier is good for people with a recessed or flat stoma
- Tap at bottom allows for easy draining of urine



M

CIN	Mfr. No.	Color	Barrier Opening	Qty.
1835594	175792	Transparent	¾ in.	5/ct
1835602	175793	Transparent	⅞ in.	5/ct
1835610	175794	Transparent	1 in.	5/ct
1835628	175795	Transparent	1⅛ in.	5/ct
1835636	175796	Transparent	1¼ in.	5/ct

HCPC Code: A4393

M May be reimbursable by Medicare Part B

Available through Cardinal Health Drop Ship program

Note: Products may not be available in all Cardinal Health Distribution Centers

## ConvaTec® - Barriers/Flanges

ConvaTec®

### Esteem Synergy® Stomahesive® Skin Barrier with Flexible Tape Collar Cut-to-Fit



- Features an innovative Adhesive Coupling Technology™
- Instead of plastic rings, a patented foam adhesive allows the pouch to be securely attached to the skin barrier on a clear plastic "Landing Zone" flange
- Stomahesive® skin barriers are designed to be gentle on peristomal skin
- Flexible tape collar

M

CIN	Mfr. No.	Description	Color	Barrier Opening	Qty.
3531092	405457	Medium	White	Up to 1⅞ in.	10/ct
3531100	405458	Large	White	Up to 2⅞ in.	10/ct

HCPC Code: A4409

ConvaTec®

### SUR-FIT Natura® Two-Piece Cut-to-Fit Stomahesive® Skin Barrier with Flange



- A comprehensive product line with the latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Stomahesive® skin barriers are ideal for people whose stoma output is semi-formed to formed
- Reliable skin barrier provides a unique formula that enables it to adhere to dry and moist skin
- Gentle on skin yet holds firmly

M

CIN	Mfr. No.	Flange Size	Size	Qty.
<b>No tape collar</b>				
2846830	401573	1¼ in.	4 x 4 in.	10/ct
2848166	401574	1½ in.	4 x 4 in.	10/ct
2847499	401575	1¾ in.	4 x 4 in.	10/ct
2848158	401576	2¼ in.	4 x 4 in.	10/ct
2846954	401577	2¾ in.	5 x 5 in.	10/ct
2847481	401906	4 in.	6 x 6 in.	5/ct

HCPC Code: A4415

ConvaTec®

### SUR-FIT Natura® Two-Piece Durahesive® Skin Barrier with CONVEX-IT® Pre-Cut



- A comprehensive product line with the latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Designed for people whose stoma output is mostly liquid
- Durahesive® skin barriers swell up or "turtleneck" to further protect the stoma and peristomal skin

M

CIN	Mfr. No.	Flange Size	Stoma Size	Qty.
<b>With flexible tape collar</b>				
3490729	413178	1¾ in.	⅝ in.	10/ct
3490737	413179	1¾ in.	¾ in.	10/ct
3490745	413180	1¾ in.	⅞ in.	10/ct
3490752	413181	1¾ in.	1 in.	10/ct
3490760	413182	1¾ in.	1⅞ in.	10/ct
3490778	413183	1¾ in.	1¼ in.	10/ct
3490786	413184	1¾ in.	1⅝ in.	10/ct
3490794	413185	2¼ in.	1½ in.	10/ct

HCPC Code: A4407



ConvaTec®

**SUR-FIT Natura® Two-Piece Durahesive® Skin Barrier Cut-to-Fit**

- A comprehensive product line with the latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Designed for people whose stoma output is mostly liquid
- Durahesive® skin barriers swell up or “turtleneck” to further protect the stoma and peristomal skin



M

CIN	Mfr. No.	Flange Size	Size	Qty.
<b>No tape collar</b>				
3490141	413154	1½ in.	4 x 4 in.	10/ct
3490182	413155	1¾ in.	4 x 4 in.	10/ct
3490190	413156	2¼ in.	4 x 4 in.	10/ct
3490216	413157	2¾ in.	5 x 5 in.	10/ct
2848364	401905	4 in.	6 x 6 in.	5/ct

HCPC Code: A4409

ConvaTec®

**SUR-FIT Natura® Two-Piece Durahesive® Skin Barrier Cut-to-Fit**

- A comprehensive product line with the latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Designed for people whose stoma output is mostly liquid



M

CIN	Mfr. No.	Flange Size	Barrier Opening	Qty.
<b>With flexible tape collar, white</b>				
3490232	413160	1½ in.	4.5 x 4.5 in.	10/ct
3490240	413161	1¾ in.	4.5 x 4.5 in.	10/ct
3490257	413162	2¼ in.	5 x 5 in.	10/ct
3490265	413163	2¾ in.	5 x 5 in.	10/ct
<b>With flexible tape collar, tan</b>				
3490281	413165	1½ in.	4 x 4 in.	10/ct
3490299	413166	1¾ in.	4 x 4 in.	10/ct
3490315	413167	2¼ in.	5 x 5 in.	10/ct
3490356	413168	2¾ in.	5 x 5 in.	10/ct

HCPC Code: A4409

ConvaTec®

**SUR-FIT Natura® Durahesive® Skin Barrier Moldable Convex**

- A comprehensive product line with the latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Designed for people whose stoma output is mostly liquid
- Durahesive® skin barriers swell up or “turtleneck” to further protect the stoma and peristomal skin
- Convex is good for people with a retracted or flush stoma and moldable for improved skin protection over traditional barriers—provides ease of use and versatility



M

CIN	Mfr. No.	Flange Size	Stoma Size	Size	Qty.
<b>With flexible tape collar, white</b>					
3462298	404592	1¾ in.	½–¾ in.	Small	10/ct
3462637	404593	1¾ in.	¾–1¼ in.	Medium	10/ct
3462645	404594	2¼ in.	1¼–1¾ in.	Large	10/ct

HCPC Code: A4407

ConvaTec®

**SUR-FIT Natura® Durahesive® Skin Barrier Moldable**

- A comprehensive product line with the latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Designed for people whose stoma output is mostly liquid
- Durahesive® skin barriers swell up or “turtleneck” to further protect the stoma and peristomal skin



M

CIN	Mfr. No.	Flange Size	Stoma Size	Size	Qty.
<b>With hydrocolloid flexible collar, white</b>					
4007522	411802	1¾ in.	¾–1¼ in.	Medium	10/ct
4021416	411804	2¼ in.	1¼–1¾ in.	Large	10/ct
4022612	411806	2¾ in.	1¾–2½ in.	X-large	10/ct

HCPC Code: A4409

ConvaTec®

**SUR-FIT Natura® Stomahesive® Skin Barrier Moldable**

- A comprehensive product line with the latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Designed for people whose stoma output is mostly liquid
- Durahesive® skin barriers swell up or “turtleneck” to further protect the stoma and peristomal skin



M

CIN	Mfr. No.	Flange Size	Stoma Size	Size	Qty.
<b>With hydrocolloid flexible collar, white</b>					
4021382	411803	1¾ in.	¾–1¼ in.	Medium	10/ct
4021390	411805	2¼ in.	1¼–1¾ in.	Large	10/ct
4022604	411807	2¾ in.	1¾–2½ in.	X-large	10/ct

HCPC Code: A4409

M May be reimbursable by Medicare Part B

Available through Cardinal Health Drop Ship program

Note: Products may not be available in all Cardinal Health Distribution Centers

## ConvaTec® - Barriers/Flanges

ConvaTec®

### SUR-FIT Natura® Two-Piece Stomahesive® Skin Barrier Pre-Cut

- A comprehensive product line with our latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Ideal for people whose stoma output is semi-formed to formed
- Reliable skin barrier provides a unique formula that enables the skin barrier to adhere to both dry and moist skin
- Gentle on skin yet holds firmly



M

CIN	Mfr. No.	Flange Size	Stoma Size	Qty.
<b>With flexible tape collar, tan</b>				
2848117	125269	1¾ in.	¾ in.	10/ct
2847903	125270	1¾ in.	7⁄8 in.	10/ct
2847721	125271	1¾ in.	1 in.	10/ct
2847648	125272	1¾ in.	1½ in.	10/ct
2847622	125273	1¾ in.	1¼ in.	10/ct
2847614	125274	1¾ in.	1⅜ in.	10/ct
2848349	125275	2¼ in.	1½ in.	10/ct
2847598	125277	2¼ in.	1¾ in.	10/ct

HCPC Code: A4414

ConvaTec®

### SUR-FIT Natura® Two-Piece Stomahesive® Skin Barrier Cut-to-Fit

- A comprehensive product line with our latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Ideal for people whose stoma output is semi-formed to formed
- Reliable skin barrier provides a unique formula that enables the skin barrier to adhere to both dry and moist skin
- Gentle on skin yet holds firmly



M

CIN	Mfr. No.	Flange Size	Size	Qty.
<b>With flexible tape collar, white</b>				
2551141	125257	1¼ in.	4 x 4 in.	10/ct
2551158	125258	1½ in.	4 x 4 in.	10/ct
2551166	125259	1¾ in.	4 x 4 in.	10/ct
2551174	125260	2¼ in.	5 x 5 in.	10/ct
2551182	125261	2¾ in.	5 x 5 in.	10/ct
<b>With flexible tape collar, tan</b>				
2551091	125262	1¼ in.	4 x 4 in.	10/ct
2551109	125263	1½ in.	4 x 4 in.	10/ct
2551117	125264	1¾ in.	4 x 4 in.	10/ct
2551125	125265	2¼ in.	5 x 5 in.	10/ct

HCPC Code: A4414

## ConvaTec® - Two-Piece Systems

ConvaTec®

### SUR-FIT Natura® Two-Piece Closed-End Pouch

- Designed for one-time usage
- No draining or uncomfortable clip to worry about
- Simply remove the pouch and place it into a disposable bag to discard
- Compatible with Stomahesive® and Durahesive® skin barriers



M

CIN	Mfr. No.	Color	Flange Size	Qty.
<b>2-sided comfort panel</b>				
3556263	413170	Opaque	1¾ in.	60/ct
3556313	413171	Opaque	2¼ in.	60/ct
3581642	413172	Opaque	2¾ in.	60/ct
2551349	401521	Opaque	1½ in.	30/ct
2551356	401522	Opaque	1¾ in.	30/ct
2551364	401523	Opaque	2¼ in.	30/ct
2551372	401524	Opaque	2¾ in.	30/ct

HCPC Code: A5054

ConvaTec®

**SUR-FIT Natura® Closed-End Pouch with Filter**

- Designed for one-time usage
- Filter is designed to reduce odor and gas buildup
- No draining or uncomfortable clip to worry about
- Simply remove the pouch and place it into a disposable bag to discard
- Compatible with Stomahesive® and Durahesive® skin barriers



M

CIN	Mfr. No.	Color	Flange Size	Qty.
<b>2-sided comfort panel</b>				
2551380	401525	Opaque	1½ in.	30/ct
2551406	401526	Opaque	1¾ in.	30/ct
2551414	401527	Opaque	2¼ in.	30/ct
2551422	401528	Opaque	2¾ in.	30/ct

HCPC Code: A4419



ConvaTec®

**SUR-FIT Natura® Two-Piece Drainable Pouch**

- A comprehensive product line with the latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Only the SUR-FIT Natura® ostomy system is available with ConvaTec® Moldable Technology™ skin barriers
- These barriers hug the individual contours of the stoma to allow for a secure, personalized fit without the need for scissors
- Pouch allows for draining



M

CIN	Mfr. No.	Flange Size	Pouch Size	Qty.
<b>1-sided comfort panel, transparent</b>				
2551307	401511	1½ in.	12 in.	10/ct
2551315	401512	1¾ in.	12 in.	10/ct
2551323	401513	2¼ in.	12 in.	10/ct
2551331	401514	2¾ in.	12 in.	10/ct
2847333	401515	4 in.	14 in.	10/ct
<b>1-sided comfort panel, opaque</b>				
2847234	401932	1¼ in.	12 in.	10/ct
2848174	401933	1½ in.	12 in.	10/ct
2847283	401934	1¾ in.	12 in.	10/ct
2848141	401935	2¼ in.	12 in.	10/ct
2847309	401936	2¾ in.	12 in.	10/ct
<b>2-sided comfort panel, opaque</b>				
2551240	401500	1¼ in.	12 in.	10/ct
2551257	401501	1½ in.	12 in.	10/ct
2551265	401502	1¾ in.	12 in.	10/ct
2551273	401503	2¼ in.	12 in.	10/ct
2551281	401504	2¾ in.	12 in.	10/ct
2551190	401505	1¼ in.	10 in.	10/ct
2551208	401506	1½ in.	10 in.	10/ct
2551216	401507	1¾ in.	10 in.	10/ct
2551224	401508	2¼ in.	10 in.	10/ct
2551232	401509	2¾ in.	10 in.	10/ct

HCPC Code: A5063

**If you don't see what you're looking for... ask us about it.**

This catalog contains just some of the items available through GoodPill. If you don't see the item you need, call **754.201.3663** and ask a team member for assistance.

M May be reimbursable by Medicare Part B

Available through Cardinal Health Drop Ship program

Note: Products may not be available in all Cardinal Health Distribution Centers

## ConvaTec® - Two-Piece Systems

### ConvaTec®

#### SUR-FIT Natura® Two-Piece Drainable Pouch with Filter

- A comprehensive product line with the latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Only the SUR-FIT Natura® ostomy system is available with ConvaTec® Moldable Technology™ skin barriers
- These barriers hug the individual contours of the stoma to allow for a secure, personalized fit without the need for scissors
- Pouch allows for draining



M

CIN	Mfr. No.	Flange Size	Pouch Size	Qty.
<b>2-sided comfort panel, opaque</b>				
3610904	411491	1¾ in.	12 in.	20/ct
3610912	411492	2¼ in.	12 in.	20/ct
3610920	411493	2¾ in.	12 in.	20/ct
3610953	411290	1¾ in.	10 in.	20/ct
3610961	411291	2¼ in.	10 in.	20/ct
3610979	411292	2¾ in.	10 in.	20/ct
<b>1-sided comfort panel, transparent</b>				
3610854	411265	1¾ in.	12 in.	10/ct
3610862	411266	2¼ in.	12 in.	10/ct
3610870	411267	2¾ in.	12 in.	10/ct

HCPC Code: A5063

### ConvaTec®

#### SUR-FIT Natura® Two-Piece Drainable Pouch with InvisiClose® with Filter

- A comprehensive product line with the latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Only the SUR-FIT Natura® ostomy system is available with ConvaTec® Moldable Technology™ skin barriers
- These barriers hug the individual contours of the stoma to allow for a secure, personalized fit without the need for scissors
- Pouch allows for draining



M

CIN	Mfr. No.	Flange Size	Pouch Size	Qty.
<b>2-sided comfort panel, opaque</b>				
3666856	411310	1¾ in.	12 in.	20/ct
3666864	411311	2¼ in.	12 in.	20/ct
3666898	411312	2¾ in.	12 in.	20 ct
<b>1-sided comfort panel, transparent</b>				
3666773	411360	1¾ in.	12 in.	10/ct
3666807	411361	2¼ in.	12 in.	10/ct

HCPC Code: A5063

### ConvaTec®

#### SUR-FIT Natura® Two-Piece Urostomy Pouch Accuseal® Tap with Valve

- A comprehensive product line with the latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Extended anti-reflux valve, comfort panels and a choice of two types of drainage taps



M

CIN	Mfr. No.	Flange Size	Size	Qty.
<b>1-sided comfort panel, transparent</b>				
2551430	401542	1¼ in.	Standard	10/ct
2551448	401543	1½ in.	Standard	10/ct
2551455	401544	1¾ in.	Standard	10/ct
2551463	401545	2¼ in.	Standard	10/ct
2551471	401546	2¾ in.	Standard	10/ct
<b>1-sided comfort panel, opaque</b>				
2551505	401549	1¾ in.	Small	10/ct
2551539	401552	1½ in.	Standard	10/ct
2551547	401553	1¾ in.	Standard	10/ct
2551554	401554	2¼ in.	Standard	10/ct

HCPC Code: A4432

### ConvaTec®

#### SUR-FIT Natura® Two-Piece Urostomy Pouch with Fold-Up Tap

- A comprehensive product line with the latest skin barrier technologies
- Plastic coupling ring allows the pouch and skin barrier to snap together for audible security
- Extended anti-reflux valve, comfort panels and a choice of two types of drainage taps



M

CIN	Mfr. No.	Flange Size	Size	Qty.
<b>1-sided comfort panel, transparent</b>				
2846780	401535	1¾ in.	Standard	10/ct
2848455	401536	2¼ in.	Standard	10/ct

HCPC Code: A4432



## Hollister - One-Piece Systems

### Hollister

#### Premier™ Drainable Pouch with Flextend® Cut-to-Fit, Flat Barrier and Lock 'n Roll® Closure

- Protects the skin over extended wearing time from irritating discharge
- Security and confidence assured with porous tape
- ComfortWear panels provide a soft, cloth-like covering between the pouch and the skin



M

CIN	Mfr. No.	Color	Barrier Opening	Pouch Size	Qty.
3613023	8531	Ultra-clear	Up to 2½ in.	12 in.	10/ct

HCPC Code: A4388

### Hollister

#### Premier™ Drainable Pouch with Flextend® Cut-to-Fit, Flat Barrier and Pouch Clamp Closure

- Protects the skin over extended wearing time from irritating discharge
- ComfortWear panels provide a soft, cloth-like covering between the pouch and the skin
- Ultimate in confidence and discreetness with odor-barrier, rustle-free film and anatomically-shaped pouches



M

CIN	Mfr. No.	Color	Barrier Opening	Pouch Size	Qty.
2316008	8631	Transparent	Up to 2½ in.	12 in.	10/ct

HCPC Code: A4388

### Hollister

#### Premier™ Drainable Pouch with Flextend® Skin Barrier Convex Cut-to-Fit and Pouch Clamp Closure

- Built-in convexity means added security and skin protection for recessed, retracted and flush stomas
- Protects the skin over extended wearing time from irritating discharge
- ComfortWear panels provide a soft, cloth-like covering between the pouch and the skin
- Ultimate in confidence and discreetness with odor-barrier, rustle-free film and anatomically-shaped pouches



M

CIN	Mfr. No.	Color	Barrier Opening	Pouch Size	Qty.
2570463	8628	Transparent	Up to 1½ in.	12 in.	5/ct

HCPC Code: A4390

### Hollister

#### Premier™ Drainable Pouch with Flextend® Skin Barrier Convex Pre-Sized and Pouch Clamp Closure

- Built-in convexity means added security and skin protection for recessed, retracted and flush stomas
- Protects the skin over extended wearing time from irritating discharge
- ComfortWear panels provide a soft, cloth-like covering between the pouch and the skin
- Ultimate in confidence and discreetness with odor-barrier, rustle-free film and anatomically-shaped pouches



M

CIN	Mfr. No.	Color	Barrier Opening	Pouch Size	Qty.
3025814	8610	Transparent	¾ in.	12 in.	5/ct
2080679	8611	Transparent	⅞ in.	12 in.	5/ct
2696359	8612	Transparent	1 in.	12 in.	5/ct
2487197	8613	Transparent	1⅛ in.	12 in.	5/ct
2360691	8614	Transparent	1¼ in.	12 in.	5/ct

HCPC Code: A4390

### Hollister

#### Premier™ Urostomy Pouch with Convex Barrier Cut-to-Fit

- Built-in convexity means added security and skin protection for recessed, retracted and flush stomas
- Protects the skin over extended wearing time from irritating discharge
- Security and confidence assured with internal non-reflux valves and a drain valve that allows regulation of flow rate
- With tape border



M

CIN	Mfr. No.	Color	Barrier Opening	Pouch Size	Qty.
2570448	8478	Transparent	Up to 1½ in.	12 in.	5/ct

HCPC Code: A4430

M May be reimbursable by Medicare Part B

Available through Cardinal Health Drop Ship program

Note: Products may not be available in all Cardinal Health Distribution Centers



## Hollister - One-Piece Systems



**Hollister**

**Premier™ Urostomy Pouch with Convex Barrier Pre-Sized**

- Built-in convexity means added security and skin protection for recessed, retracted and flush stomas
- Protects the skin over extended wearing time from irritating discharge
- Security and confidence assured with internal non-reflux valves and a drain valve that allows regulation of flow rate

**M**

CIN	Mfr. No.	Color	Stoma Size	Pouch Size	Qty.
2193696	8481	Transparent	5/8 in.	9 in.	5/ct
2590263	8482	Transparent	3/4 in.	9 in.	5/ct
2294965	8483	Transparent	7/8 in.	9 in.	5/ct
2561884	8484	Transparent	1 in.	9 in.	5/ct
2181915	8485	Transparent	1 1/8 in.	9 in.	5/ct
2431864	8486	Transparent	1 1/4 in.	9 in.	5/ct
2570455	8487	Transparent	1 3/8 in.	9 in.	5/ct
3233343	8488	Transparent	1 1/2 in.	9 in.	5/ct
2982106	8489	Transparent	1 3/4 in.	9 in.	5/ct

**HCPC Code: A4430**



## Hollister - Barriers/Flanges

**Hollister**

**CenterPointLock™ FlexWear® Flat Skin Barrier with Floating Flange and Tape Border Filter**

- Protects skin from corrosive drainage with the standard wear FlexWear® skin barrier
- Eliminates pressure on tender post-operative abdomens with the exclusive Floating Flange
- Confidence and security with the unique CenterPointLock™ feature—prevents unintentional pouch/skin barrier separation



**M**


CIN	Mfr. No.	Flange Size	Barrier Opening	Qty.
3041498	3727	1 1/2 in.	Up to 1 in.	5/ct
1948322	3722	1 3/4 in.	Up to 1 1/4 in.	5/ct
2539526	3723	2 1/4 in.	Up to 1 3/4 in.	5/ct
2288884	3724	2 3/4 in.	Up to 2 1/4 in.	5/ct
2342012	3726	4 in.	Up to 3 1/2 in.	5/ct

**HCPC Code: A4414**

**Hollister**

**CenterPointLock™ Convex Skin Barrier with Floating Flange**

- Built-in convexity means added security and skin protection for recessed, retracted and flush stomas—cut-to-fit openings
- Protects the skin from irritating discharge over extended wearing time
- Security assured with patented CenterPointLock™ feature—prevents accidental separation of pouch from skin barrier
- With tape border, porous cloth adhesive



**M**

CIN	Mfr. No.	Flange Size	Barrier Opening	Qty.
2538353	8778	2 1/4 in.	Up to 1 1/2 in.	5/ct

**HCPC Code: A4407**



Hollister - Barriers/Flanges



Hollister

**CenterPointLock™ Flextend® Skin Barrier Cut-to-Fit Flat with Floating Flange**

- Protects the skin over extended wearing time from irritating discharge
- Confidence and security with the unique CenterPointLock™ feature—prevents unintentional pouch/skin barrier separation
- Eliminates pressure on tender post-operative abdomens with the exclusive Floating Flange
- With tape border, porous cloth adhesive

M

CIN	Mfr. No.	Flange Size	Barrier Opening	Qty.
2599066	8722	1¾ in.	Up to 1¼ in.	5/ct
2256741	8723	2¼ in.	Up to 1¾ in.	5/ct
2531283	8724	2¾ in.	Up to 2¼ in.	5/ct

HCPC Code: A4407



Hollister

**CenterPointLock™ Hollihesive Skin Barrier Cut-to-Fit with Floating Flange**

- Eliminates pressure on tender post-operative abdomens with the exclusive Floating Flange
- Confidence and security with the unique CenterPointLock™ feature—prevents unintentional pouch/skin barrier separation
- With tape border, porous paper adhesive

M

CIN	Mfr. No.	Flange Size	Barrier Opening	Qty.
1919075	3702	1¾ in.	Up to 1 in.	5/ct
1948462	3703	2¼ in.	Up to 1½ in.	5/ct
1919091	3704	2¾ in.	Up to 2 in.	5/ct
1919109	3706	4 in.	Up to 3¾ in.	5/ct

HCPC Code: A4414, A4415

Hollister

**New Image® Flextend® Convex Skin Barrier Cut-to-Fit with Floating Flange**




- Built-in convexity means added security and skin protection for recessed, retracted and flush stomas—cut-to-fit openings
- Eliminates pressure on tender post-operative abdomens with the exclusive Floating Flange
- Protects the skin over extended wearing time from irritating discharge

M

CIN	Mfr. No.	Flange Size	Barrier Opening	Qty.
<b>With tape border</b>				
3674769	14802	1¾ in.	Up to 1 in.	5/ct
3638343	14803	2¼ in.	Up to 1½ in.	5/ct
3651346	14804	2¾ in.	Up to 2 in.	5/ct

HCPC Code: A4407



Hollister

**New Image® Flextend® Convex Skin Barrier Pre-Sized with Floating Flange**

- Built-in convexity means added security and skin protection for recessed, retracted and flush stomas—pre-cut openings
- Eliminates pressure on tender post-operative abdomens with the exclusive Floating Flange
- Protects the skin over extended wearing time from irritating discharge
- Flextend skin barrier is the barrier of choice for urostomies and ileostomies

M

CIN	Mfr. No.	Flange Size	Barrier Opening	Qty.
<b>With tape border</b>				
3571825	14903	1¾ in.	¾ in.	5/ct
3586393	14904	1¾ in.	1 in.	5/ct
3895471	14906	2¼ in.	1¼ in.	5/ct

HCPC Code: A4407

M May be reimbursable by Medicare Part B

Available through Cardinal Health Drop Ship program

Note: Products may not be available in all Cardinal Health Distribution Centers

## Hollister - Barriers/Flanges

### Hollister

#### New Image® Flexlend® Skin Barrier Cut-to-Fit with Floating Flange



- Confidence assured with the **Floating Flange**—provides utmost security with the patented locking system that audibly “clicks” into place
- Eliminates pressure on tender post-operative abdomens with the exclusive **Floating Flange**
- Protects the skin over extended wearing time
- Flexlend® skin barrier with tapered edges eliminates the need for tape

M

CIN	Mfr. No.	Flange Size	Barrier Opening	Qty.
<b>Without tape border</b>				
3146479	15602	1¾ in.	Up to 1¼ in.	5/ct
3037546	15603	2¼ in.	Up to 1¾ in.	5/ct
3149234	15604	2¾ in.	Up to 2¼ in.	5/ct

HCPC Code: A4414

### Hollister

#### New Image® Flexlend® Skin Barrier Cut-to-Fit Flat with Integrated Floating Flange



- Confidence assured with the **Floating Flange**—provides utmost security with the patented locking system that audibly “clicks” into place
- Protects the skin over extended wearing time from irritating discharge
- Eliminates pressure on tender post-operative abdomens with the exclusive **Floating Flange**

M

CIN	Mfr. No.	Flange Size	Barrier Opening	Qty.
<b>With tape border</b>				
3461894	14602	1¾ in.	Up to 1¼ in.	5/ct
3461902	14603	2¼ in.	Up to 1¾ in.	5/ct
3461928	14604	2¾ in.	Up to 2¼ in.	5/ct

HCPC Code: A4414

### Hollister

#### New Image® FlexWear® Skin Barrier Cut-to-Fit with Floating Flange



- Confidence assured with the **Floating Flange**—provides utmost security with the patented locking system that audibly “clicks” into place
- Protects the skin over extended wearing time from irritating discharge
- Eliminates pressure on tender post-operative abdomens with the exclusive **Floating Flange**

M

CIN	Mfr. No.	Flange Size	Barrier Opening	Qty.
<b>With tape border</b>				
3461837	14202	1¾ in.	Up to 1¼ in.	5/ct
3461852	14203	2¼ in.	Up to 1¾ in.	5/ct
<b>Without tape border</b>				
3136900	15202	1¾ in.	Up to 1¼ in.	5/ct
3037553	15203	2¼ in.	Up to 1¾ in.	5/ct
3036530	15204	2¾ in.	Up to 2¼ in.	5/ct

HCPC Code: A4414, A4415

## Hollister - Two-Piece Systems

### Hollister

#### CenterPointLock™ Closed Pouch with Filter

- Security assured with patented CenterPointLock™ feature—prevents accidental separation of pouch from skin barrier
- Four layer odor-barrier, rustle-free film is soft and silent
- Security and confidence achieved by highly effective deodorizing and venting filter made of activated carbon granules



M

CIN	Mfr. No.	Color	Flange Size	Pouch Size	Qty.
2539534	3343	Beige	2¼ in.	9 in.	15/ct
2538296	3344	Beige	2¾ in.	9 in.	15/ct

HCPC Code: A4419

### Hollister

#### CenterPointLock™ Urostomy Pouch

- Security and confidence assured with internal non-reflux valves and a drain valve that allows regulation of flow rate
- Confidence and security with the unique CenterPointLock™ feature—prevents unintentional pouch/skin barrier separation
- Ultimate in confidence with odor-barrier, rustle-free film—four layer odor-barrier, rustle-free film is soft and silent
- Includes belt tabs, anti-reflux valve and adjustable drain valve



M

CIN	Mfr. No.	Flange Size	Pouch Size	Qty.
1919257	3902	1¾ in.	9 in.	10/ct
1919265	3903	2¼ in.	9 in.	10/ct
1919273	3904	2¾ in.	9 in.	10/ct

HCPC Code: A4432

### Hollister

#### CenterPointLock™ Drainable Pouch with Pouch Clamp Closure Filter

- Security assured with patented CenterPointLock™ feature—prevents accidental separation of pouch from skin barrier
- Four layer odor-barrier, rustle-free film is soft and silent
- Peace of mind knowing that the pouch is discreet and made of the highest quality materials



M

CIN	Mfr. No.	Color	Flange Size	Pouch Size	Qty.
1919174	3807	Transparent	1½ in.	11 in.	10/ct
1919133	3802	Transparent	1¾ in.	11 in.	10/ct
1919141	3803	Transparent	2¼ in.	11 in.	10/ct
1919158	3804	Transparent	2¾ in.	11 in.	10/ct
1919166	3806	Transparent	4 in.	11 in.	10/ct
1919182	3812	Beige	1¾ in.	11 in.	10/ct
1919190	3813	Beige	2¼ in.	11 in.	10/ct
1919208	3814	Beige	2¾ in.	11 in.	10/ct

HCPC Code: A5063

### Hollister

#### New Image® Drainable Pouch with Pouch Clamp Closure

- ComfortWear panels provide a soft, cloth-like covering between the pouch and the skin
- Confidence and security with the unique CenterPointLock™ feature—prevents unintentional pouch/skin barrier separation
- Provides maximum protection, utilizes an odor-barrier material and are rustle-free



M

CIN	Mfr. No.	Flange Size	Pouch Size	Qty.
<b>Panel on body side only, transparent</b>				
3146461	18102	1¾ in.	12 in.	10/ct
3105822	18103	2¼ in.	12 in.	10/ct
3036563	18104	2¾ in.	12 in.	10/ct
<b>Panel on both sides, beige</b>				
3337979	18122	1¾ in.	12 in.	10/ct
3260353	18123	2¼ in.	12 in.	10/ct
3149242	18124	2¾ in.	12 in.	10/ct

HCPC Code: A5063

M May be reimbursable by Medicare Part B

Available through Cardinal Health Drop Ship program

Note: Products may not be available in all Cardinal Health Distribution Centers

## Hollister - Two-Piece Systems

### Hollister

#### New Image® Drainable Pouch with Lock 'n Roll® Closure without Filter

- Easy to close—simply roll it three times in the same direction, press from the middle out and **feel the seal**
- Easy to empty—lightly pinch the sides of the leak-proof outlet for easy control of content flow; firm reinforcing strips shut as soon as you stop pinching
- Easy to clean—simply clean the inside of the pouch tail by wiping it with a damp tissue



M

CIN	Mfr. No.	Color	Flange Size	Pouch Size	Qty.
<b>Panel on body side only</b>					
3571833	18132	Transparent	1¾ in.	12 in.	10/ct
3571841	18133	Transparent	2¼ in.	12 in.	10/ct
3571858	18134	Transparent	2¾ in.	12 in.	10/ct

HCPC Code: A5063



### Hollister

#### New Image® Drainable Pouch with Lock 'n Roll® Closure

- Easy to close—simply roll it three times in the same direction, press from the middle out and **feel the seal**
- Easy to empty—lightly pinch the sides of the leak-proof outlet for easy control of content flow; firm reinforcing strips shut as soon as you stop pinching
- Easy to clean—simply clean the inside of the pouch tail by wiping it with a damp tissue



M

CIN	Mfr. No.	Flange Size	Pouch Size	Qty.
<b>Panel on body side only, transparent</b>				
3525441	18192	1¾ in.	12 in.	10/ct
3525458	18193	2¼ in.	12 in.	10/ct
3525466	18194	2¾ in.	12 in.	10/ct
<b>Panel on both sides, beige</b>				
3551777	18182	1¾ in.	12 in.	10/ct
3551793	18183	2¼ in.	12 in.	10/ct
3548021	18184	2¾ in.	12 in.	10/ct

HCPC Code: A4425

### Hollister

#### New Image® Urostomy Pouch

- Security and confidence assured with internal non-reflux valves and a drain valve that allows regulation of flow rate
- Ultimate in comfort with soft ComfortWear panels against the skin
- Ultimate in confidence and discreetness with odor-barrier, rustle-free film and anatomically-shaped pouches
- Includes belt tabs, anti-reflux valve and adjustable drain valve



M

CIN	Mfr. No.	Flange Size	Pouch Size	Qty.
<b>Panel on body side only, transparent</b>				
2982130	18402	1¾ in.	9 in.	10/ct
3302627	18403	2¼ in.	9 in.	10/ct
3155504	18404	2¾ in.	9 in.	10/ct

HCPC Code: A4432